

Product datasheet for **TP720999XL**

SF20 (MYDGF) (NM_019107) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human chromosome 19 open reading frame 10 (C19orf10)
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Ser33-Leu173
Tag:	N-His
Predicted MW:	18 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Lyophilized from a 0.2 um filtered solution of 20mMTris-HCl, 150mMNaCl,pH8.0.
Endotoxin:	Endotoxin level is < 0.1 ng/μg of protein (< 1 EU/μg)
Reconstitution Method:	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in ddH ₂ O. It is not recommended to reconstitute a concentration less than 100 μg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Storage:	Lyophilized protein should be stored at < -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_061980
Locus ID:	56005
UniProt ID:	Q969H8
RefSeq Size:	1067
Cytogenetics:	19p13.3
RefSeq ORF:	519
Synonyms:	C19orf10; EUROIMAGE1875335; IL25; IL27; IL27w; R33729_1; SF20



[View online »](#)

Summary: The protein encoded by this gene was previously thought to support proliferation of lymphoid cells and was considered an interleukin. However, this activity has not been reproducible and the function of this protein is currently unknown. [provided by RefSeq, Jul 2008]

Protein Families: Secreted Protein