

Product datasheet for **TP721004XL**

ZNHIT1 (NM_006349) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human zinc finger, HIT-type containing 1 (ZNHIT1)
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Met1-Val154
Tag:	N-His
Predicted MW:	19.6 kDa
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 50% Glycerol, pH 7.4.
Endotoxin:	Endotoxin level is < 0.1 ng/μg of protein (< 1 EU/μg)
Storage:	Store at < -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_006340
Locus ID:	10467
UniProt ID:	O43257
RefSeq Size:	1208
Cytogenetics:	7q22.1
RefSeq ORF:	462
Synonyms:	CG1I; ZNFN4A1
Summary:	Seems to play a role in p53-mediated apoptosis induction. Binds to NR1D2 and relieves it of its inhibitory effect on the transcription of APOC3 without affecting its DNA-binding activity. [UniProtKB/Swiss-Prot Function]



[View online »](#)