

Product datasheet for **TP721007XL**

UPRT (NM_145052) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human uracil phosphoribosyltransferase (FUR1) homolog (S. cerevisiae) (UPRT), transcript variant 1
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Met1-Asp309
Tag:	N-His
Predicted MW:	35.9 kDa
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_659489
Locus ID:	139596
UniProt ID:	Q96BW1 , A8KAF9
RefSeq Size:	2512
Cytogenetics:	Xq13.3
RefSeq ORF:	927
Synonyms:	FUR1; UPP
Summary:	This gene encodes uracil phosphoribosyltransferase, which catalyzes the conversion of uracil and 5-phosphoribosyl-1-R-diphosphate to uridine monophosphate (UMP). This reaction is an important part of nucleotide metabolism, specifically the pyrimidine salvage pathway. The enzyme localizes to the nucleus and cytoplasm. The protein is a potential target for rational design of drugs to treat parasitic infections and cancer. [provided by RefSeq, Nov 2009]



[View online »](#)

Protein Families: Druggable Genome

Protein Pathways: Metabolic pathways, Pyrimidine metabolism