

# **Product datasheet for TP721050XL**

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Flt3 ligand (FLT3LG) (NM\_001459) Human Recombinant Protein

#### **Product data:**

**Product Type:** Recombinant Proteins

**Description:** Purified recombinant protein of Human fms-related tyrosine kinase 3 ligand (FLT3LG),

transcript variant 3

Species: Human Expression Host: HEK293

Expression cDNA Clone

or AA Sequence:

Thr27-Pro184

Tag: C-His

Predicted MW: 18.9 kDa

Concentration: lot specific

**Purity:** >95% as determined by SDS-PAGE and Coomassie blue staining

Buffer: Provided lyophilized from a 0.2 μm filtered solution of 20 mM Tris-HCl, 150 mM NaCl

Endotoxin: Endotoxin level is < 0.1 ng/μg of protein (< 1 EU/μg)

**Reconstitution Method:** Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the

lyophilized protein in ddH2O. It is not recommended to reconstitute a concentration less than 100  $\mu$ g/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Storage: Store at -80°C.

**Stability:** Stable for at least 6 months from date of receipt under proper storage and handling

conditions.

**RefSeq:** NP 001450

 Locus ID:
 2323

 UniProt ID:
 P49771

 RefSeq Size:
 1074

Cytogenetics: 19q13.33

RefSeq ORF: 705

Synonyms: FL; FLG3L; FLT3L





## Flt3 ligand (FLT3LG) (NM\_001459) Human Recombinant Protein - TP721050XL

Summary: Dendritic cells (DCs) provide the key link between innate and adaptive immunity by

recognizing pathogens and priming pathogen-specific immune responses. FLT3LG controls the development of DCs and is particularly important for plasmacytoid DCs and CD8 (see MIM

186910)-positive classical DCs and their CD103 (ITGAE; MIM 604682)-positive tissue

counterparts (summary by Sathaliyawala et al., 2010 [PubMed 20933441]).[supplied by OMIM,

Jan 2011]

**Protein Families:** Druggable Genome, ES Cell Differentiation/IPS, Secreted Protein, Transmembrane

Protein Pathways: Cytokine-cytokine receptor interaction, Hematopoietic cell lineage, Pathways in cancer