

## Product datasheet for **TP721063L**

### IL18R Beta (IL18RAP) (NM\_003853) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human interleukin 18 receptor accessory protein (IL18RAP)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Phe20-Arg356
Tag:	C-His
Predicted MW:	39.4 kDa
Concentration:	lot specific
Purity:	Greater than 95% as determined by RP-HPLC and reducing SDS-PAGE.
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)
Reconstitution Method:	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in ddH <sub>2</sub> O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Storage:	Store at -80°C.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	<a href="#">NP_003844</a>
Locus ID:	8807
UniProt ID:	<a href="#">O95256</a>
RefSeq Size:	2672
Cytogenetics:	2q12.1
RefSeq ORF:	1797
Synonyms:	ACPL; CD218b; CDw218b; IL-1R-7; IL-1R7; IL-1RAcPL; IL-18R-beta; IL-18RAcP; IL-18Rbeta; IL18RB



[View online »](#)

**Summary:**

The protein encoded by this gene is an accessory subunit of the heterodimeric receptor for interleukin 18 (IL18), a proinflammatory cytokine involved in inducing cell-mediated immunity. This protein enhances the IL18-binding activity of the IL18 receptor and plays a role in signaling by IL18. Mutations in this gene are associated with Crohn's disease and inflammatory bowel disease, and susceptibility to celiac disease and leprosy. Alternatively spliced transcript variants of this gene have been described, but their full-length nature is not known. [provided by RefSeq, Feb 2014]

**Protein Families:**

Druggable Genome, Transmembrane

**Protein Pathways:**

Cytokine-cytokine receptor interaction