

Product datasheet for **TP721085**

Aldolase C (ALDOC) (NM_005165) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human aldolase C, fructose-bisphosphate (ALDOC)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Phe2-Tyr364
Tag:	C-His
Predicted MW:	40.3 kDa
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 μ m filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/ μ g of protein (< 1 EU/ μ g)
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_005156
Locus ID:	230
UniProt ID:	P09972 , A0A024QZ64
RefSeq Size:	1665
Cytogenetics:	17q11.2
RefSeq ORF:	1092
Synonyms:	ALDC
Summary:	This gene encodes a member of the class I fructose-biphosphate aldolase gene family. Expressed specifically in the hippocampus and Purkinje cells of the brain, the encoded protein is a glycolytic enzyme that catalyzes the reversible aldol cleavage of fructose-1,6-biphosphate and fructose 1-phosphate to dihydroxyacetone phosphate and either glyceraldehyde-3-phosphate or glyceraldehyde, respectively. [provided by RefSeq, Jul 2008]



[View online »](#)

Protein Pathways: Fructose and mannose metabolism, Glycolysis / Gluconeogenesis, Metabolic pathways, Pentose phosphate pathway