

Product datasheet for **TP721090XL**

ADAMDEC1 (NM_014479) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human ADAM-like, decysin 1 (ADAMDEC1), transcript variant 1
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Lys206-Glu470
Tag:	C-His
Predicted MW:	50.4 kDa
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_055294
Locus ID:	27299
UniProt ID:	O15204 , B7Z6V5
RefSeq Size:	2348
Cytogenetics:	8p21.2
RefSeq ORF:	1410
Synonyms:	M12.219
Summary:	This encoded protein is thought to be a secreted protein belonging to the disintegrin metalloproteinase family. Its expression is upregulated during dendritic cells maturation. This protein may play an important role in dendritic cell function and their interactions with germinal center T cells. [provided by RefSeq, Jul 2008]
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)