

Product datasheet for **TP721098XL**

PPT1 (NM_000310) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human palmitoyl-protein thioesterase 1 (PPT1), transcript variant 1
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Asp28-Gly306
Tag:	C-His
Predicted MW:	32.3 kDa
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_000301
Locus ID:	5538
UniProt ID:	P50897
RefSeq Size:	2504
Cytogenetics:	1p34.2
RefSeq ORF:	918
Synonyms:	CLN1; INCL; PPT



[View online »](#)

Summary:	The protein encoded by this gene is a small glycoprotein involved in the catabolism of lipid-modified proteins during lysosomal degradation. The encoded enzyme removes thioester-linked fatty acyl groups such as palmitate from cysteine residues. Defects in this gene are a cause of infantile neuronal ceroid lipofuscinosis 1 (CLN1, or INCL) and neuronal ceroid lipofuscinosis 4 (CLN4). Two transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Dec 2008]
Protein Families:	Druggable Genome
Protein Pathways:	Fatty acid elongation in mitochondria, Lysosome, Metabolic pathways