

Product datasheet for **TP721099XL**

CPVL (NM_031311) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human carboxypeptidase, vitellogenic-like (CPVL), transcript variant 1
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Leu23-Gly476
Tag:	C-His
Predicted MW:	52.9 kDa
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_112601
Locus ID:	54504
UniProt ID:	Q9H3G5 , A0A024RA40
RefSeq Size:	1691
Cytogenetics:	7p14.3
RefSeq ORF:	1428
Synonyms:	HVLP
Summary:	The protein encoded by this gene is a carboxypeptidase and bears strong sequence similarity to serine carboxypeptidases. Carboxypeptidases are a large class of proteases that act to cleave a single amino acid from the carboxy termini of proteins or peptides. The exact function of this protein, however, has not been determined. [provided by RefSeq, Jan 2017]
Protein Families:	Druggable Genome, Protease



[View online »](#)