

Product datasheet for TP721104L

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Apolipoprotein A II (APOA2) (NM_001643) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Purified recombinant protein of Human apolipoprotein A-II (APOA2)

Species: Human Expression Host: HEK293

Expression cDNA Clone

Gln24-Gln100

or AA Sequence:

Tag: C-His

Predicted MW: 9.74 kDa

Concentration: lot specific

Purity: >95% as determined by SDS-PAGE and Coomassie blue staining

Buffer: Provided lyophilized from a 0.2 μm filtered solution of 20 mM Tris-HCl, 150 mM NaCl

Bioactivity: Surface Plasmon Ressonance (SPR) (PMID: 27514935)

Endotoxin: Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)

Reconstitution Method: Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the

lyophilized protein in ddH2O. It is not recommended to reconstitute a concentration less than 100 μ g/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Storage: Store at -80°C.

Stability: Stable for at least 6 months from date of receipt under proper storage and handling

conditions.

RefSeq: NP 001634

Locus ID: 336

 UniProt ID:
 P02652

RefSeq Size: 473

Cytogenetics: 1q23.3

RefSeq ORF: 300

Synonyms: Apo-All; ApoA-Il; apoAll





Apolipoprotein A II (APOA2) (NM_001643) Human Recombinant Protein - TP721104L

Summary: This gene encodes apolipoprotein (apo-) A-II, which is the second most abundant protein of

the high density lipoprotein particles. The protein is found in plasma as a monomer, homodimer, or heterodimer with apolipoprotein D. Defects in this gene may result in apolipoprotein A-II deficiency or hypercholesterolemia. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome, Secreted Protein

Protein Pathways: PPAR signaling pathway