

Product datasheet for **TP721113XL**

IFNW1 (NM_002177) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human interferon, omega 1 (IFNW1)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Leu22-Ser195
Tag:	C-His
Predicted MW:	21.1 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)
Reconstitution Method:	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in ddH ₂ O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Storage:	Store at -20°C.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_002168
Locus ID:	3467
UniProt ID:	P05000
RefSeq Size:	1514
Cytogenetics:	9p21.3
RefSeq ORF:	585



[View online »](#)

- Summary:** The protein encoded by this gene is an interferon and possesses antiviral activity. The encoded protein binds to the interferon alpha/beta receptor but not to the interferon gamma receptor. This intronless gene has several pseudogenes spread throughout the genome. [provided by RefSeq, Nov 2015]
- Protein Families:** Druggable Genome, Secreted Protein
- Protein Pathways:** Cytokine-cytokine receptor interaction, Jak-STAT signaling pathway, RIG-I-like receptor signaling pathway