

Product datasheet for **TP721114**

CEACAM21 (NM_001098506) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human carcinoembryonic antigen-related cell adhesion molecule 21 (CEACAM21), transcript variant 1
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Trp35-Gly240
Tag:	C-His
Predicted MW:	24.2 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)
Reconstitution Method:	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in ddH ₂ O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Storage:	Store at -80°C.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_001091976
Locus ID:	90273
UniProt ID:	Q3KPI0 , A0A0G2JSC8
RefSeq Size:	1378
Cytogenetics:	19q13.2
RefSeq ORF:	879
Synonyms:	CEACAM3; R29124_1
Protein Families:	Transmembrane



[View online »](#)