

Product datasheet for **TP721116XL**

BTN3A2 (NM_007047) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human butyrophilin, subfamily 3, member A2 (BTN3A2), transcript variant 1
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Gln30-Trp248
Tag:	C-His
Predicted MW:	24.6 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)
Reconstitution Method:	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in ddH ₂ O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Storage:	Store at -80°C.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_008978
Locus ID:	11118
UniProt ID:	P78410 , A8K4B5
RefSeq Size:	3776
Cytogenetics:	6p22.2
RefSeq ORF:	1005
Synonyms:	BT3.2; BTF4; BTN3.2; CD277



[View online »](#)

Summary:

This gene encodes a member of the immunoglobulin superfamily, which resides in the juxtatelomeric region of the major histocompatibility class 1 locus and is clustered with the other family members on chromosome 6. The encoded protein may be involved in the adaptive immune response. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jun 2013]

Protein Families:

Druggable Genome, Transmembrane