

## Product datasheet for **TP721117**

### ERp57 (PDIA3) (NM\_005313) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human protein disulfide isomerase family A, member 3 (PDIA3)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Ser25-Leu505
Tag:	C-His
Predicted MW:	55.3 kDa
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	<a href="#">NP_005304</a>
Locus ID:	2923
UniProt ID:	<a href="#">P30101</a> , <a href="#">V9HVV3</a>
RefSeq Size:	3060
Cytogenetics:	15q15.3
RefSeq ORF:	1515
Synonyms:	ER60; ERp57; ERp60; ERp61; GRP57; GRP58; HEL-S-93n; HEL-S-269; HsT17083; P58; PI-PLC



[View online »](#)

**Summary:**

This gene encodes a protein of the endoplasmic reticulum that interacts with lectin chaperones calreticulin and calnexin to modulate folding of newly synthesized glycoproteins. The protein was once thought to be a phospholipase; however, it has been demonstrated that the protein actually has protein disulfide isomerase activity. It is thought that complexes of lectins and this protein mediate protein folding by promoting formation of disulfide bonds in their glycoprotein substrates. This protein also functions as a molecular chaperone that prevents the formation of protein aggregates. [provided by RefSeq, Dec 2016]

**Protein Families:**

Druggable Genome

**Protein Pathways:**

Antigen processing and presentation