

## Product datasheet for **TP721119L**

### **FSTL1 (NM\_007085) Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Purified recombinant protein of Human follistatin-like 1 (FSTL1)
<b>Species:</b>	Human
<b>Expression Host:</b>	HEK293
<b>Expression cDNA Clone or AA Sequence:</b>	Glu21-Ile308
<b>Tag:</b>	C-His
<b>Predicted MW:</b>	33.5 kDa
<b>Concentration:</b>	lot specific
<b>Purity:</b>	>95% as determined by SDS-PAGE and Coomassie blue staining
<b>Buffer:</b>	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
<b>Endotoxin:</b>	Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)
<b>Reconstitution Method:</b>	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in ddH <sub>2</sub> O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
<b>RefSeq:</b>	<a href="#">NP_009016</a>
<b>Locus ID:</b>	11167
<b>UniProt ID:</b>	<a href="#">Q12841</a>
<b>RefSeq Size:</b>	3840
<b>Cytogenetics:</b>	3q13.33
<b>RefSeq ORF:</b>	924
<b>Synonyms:</b>	FRP; FSL1; MIR198; OCC-1; OCC1; tsc36



[View online »](#)

**Summary:** This gene encodes a protein with similarity to follistatin, an activin-binding protein. It contains an FS module, a follistatin-like sequence containing 10 conserved cysteine residues. This gene product is thought to be an autoantigen associated with rheumatoid arthritis. [provided by RefSeq, Jul 2008]

**Protein Families:** Secreted Protein