

## Product datasheet for TP721150M

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## CGRP (CALCB) (NM 000728) Human Recombinant Protein

**Product data:** 

**Product Type:** Recombinant Proteins

**Description:** Purified recombinant protein of Human calcitonin-related polypeptide beta (CALCB)

Species: Human
Expression Host: E. coli

**Expression cDNA Clone** 

Met1-Ala127

or AA Sequence:

Tag: N-GST

Predicted MW: 40.9 kDa

**Purity:** >95% as determined by SDS-PAGE and Coomassie blue staining

Buffer: Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl

Endotoxin: Endotoxin level is < 0.1 ng/μg of protein (< 1 EU/μg)

Storage: Store at -80°C.

**Stability:** Stable for at least 3 months from date of receipt under proper storage and handling

conditions.

**RefSeq:** NP 000719

Locus ID: 797

 UniProt ID:
 P10092

 RefSeq Size:
 1031

 Cytogenetics:
 11p15.2

RefSeq ORF: 381

**Synonyms:** CALC2; CGRP-II; CGRP2

Summary: CGRP induces vasodilation. It dilates a variety of vessels including the coronary, cerebral and

systemic vasculature. Its abundance in the CNS also points toward a neurotransmitter or

neuromodulator role.[UniProtKB/Swiss-Prot Function]

**Protein Families:** Druggable Genome, Secreted Protein

