

## Product datasheet for **TP721157XL**

### IL13 (NM\_008355) Mouse Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse interleukin 13 (IL13)
Species:	Mouse
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Ser26-Phe131
Tag:	Tag Free
Predicted MW:	11.7 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)
Reconstitution Method:	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in ddH <sub>2</sub> O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Storage:	Store at -80°C.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	<a href="#">NP_032381</a>
Locus ID:	16163
UniProt ID:	<a href="#">P20109</a>
RefSeq Size:	1217
Cytogenetics:	11 31.98 cM
RefSeq ORF:	393
Synonyms:	IL-13



[View online »](#)

**Summary:**

Cytokine. Inhibits inflammatory cytokine production. Synergizes with IL2 in regulating interferon-gamma synthesis. May be critical in regulating inflammatory and immune responses (By similarity). Positively regulates IL31RA expression in macrophages (PubMed:25847241).[UniProtKB/Swiss-Prot Function]