

## Product datasheet for **TP721163M**

### NGAL (LCN2) (NM\_005564) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human lipocalin 2 (LCN2)
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Gln21-Gly198
Tag:	C-His
Predicted MW:	21.7 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)
Storage:	Store at -80°C.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	<a href="#">NP_005555</a>
Locus ID:	3934
UniProt ID:	<a href="#">P80188</a>
RefSeq Size:	840
Cytogenetics:	9q34.11
RefSeq ORF:	594
Synonyms:	24p3; MSFI; NGAL; p25



[View online »](#)

**Summary:**

This gene encodes a protein that belongs to the lipocalin family. Members of this family transport small hydrophobic molecules such as lipids, steroid hormones and retinoids. The protein encoded by this gene is a neutrophil gelatinase-associated lipocalin and plays a role in innate immunity by limiting bacterial growth as a result of sequestering iron-containing siderophores. The presence of this protein in blood and urine is an early biomarker of acute kidney injury. This protein is thought to be involved in multiple cellular processes, including maintenance of skin homeostasis, and suppression of invasiveness and metastasis. Mice lacking this gene are more susceptible to bacterial infection than wild type mice. [provided by RefSeq, Sep 2015]

**Protein Families:**

Secreted Protein