

Product datasheet for **TP721164XL**

GABARAP (NM_007278) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human GABA(A) receptor-associated protein (GABARAP)
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Met1-Gln116
Tag:	N-His&C-Fc
Predicted MW:	43.2 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)
Storage:	Store at -80°C.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_009209
Locus ID:	11337
UniProt ID:	O95166 , Q6IAW1
RefSeq Size:	924
Cytogenetics:	17p13.1
RefSeq ORF:	351
Synonyms:	ATG8A; GABARAP-a; MM46



[View online »](#)

Summary: Gamma-aminobutyric acid A receptors [GABA(A) receptors] are ligand-gated chloride channels that mediate inhibitory neurotransmission. This gene encodes GABA(A) receptor-associated protein, which is highly positively charged in its N-terminus and shares sequence similarity with light chain-3 of microtubule-associated proteins 1A and 1B. This protein clusters neurotransmitter receptors by mediating interaction with the cytoskeleton. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome

Protein Pathways: Regulation of autophagy