

Product datasheet for **TP721166**

ZBTB9 (NM_152735) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human zinc finger and BTB domain containing 9 (ZBTB9)
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Met1-Lys473
Tag:	C-His
Predicted MW:	51.7 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)
Reconstitution Method:	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in ddH ₂ O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Storage:	Store at -80°C.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_689948
Locus ID:	221504
UniProt ID:	Q96C00 , A0A1U9X8U5
RefSeq Size:	2730
Cytogenetics:	6p21.32
RefSeq ORF:	1419
Synonyms:	ZNF919
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]



[View online »](#)

Protein Families: Transcription Factors