

# Product datasheet for TP721170M

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Fibronectin (FN1) (NM 054034) Human Recombinant Protein

#### **Product data:**

**Product Type: Recombinant Proteins** 

Description: Purified recombinant protein of Human fibronectin 1 (FN1), transcript variant 7

Species: Human E. coli **Expression Host:** 

**Expression cDNA Clone** 

or AA Sequence:

Pro1270-Ser1546&Ala1721-Thr2016

Tag: Tag Free **Predicted MW:** 62.7 kDa **Concentration:** 

**Purity:** >95% as determined by SDS-PAGE and Coomassie blue staining

**Buffer:** Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl

**Endotoxin:** Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)

**Reconstitution Method:** Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the

> lyophilized protein in ddH2O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Store at -80°C. Storage:

Stability: Stable for at least 6 months from date of receipt under proper storage and handling

conditions.

lot specific

NP 473375 RefSeq:

Locus ID: 2335 **UniProt ID:** P02751 RefSeg Size: 2402 Cytogenetics: 2q35 RefSeq ORF: 1971

CIG; ED-B; FINC; FN; FNZ; GFND; GFND2; LETS; MSF; SMDCF Synonyms:





#### Fibronectin (FN1) (NM\_054034) Human Recombinant Protein - TP721170M

**Summary:** This gene encodes fibronectin, a glycoprotein present in a soluble dimeric form in plasma,

> and in a dimeric or multimeric form at the cell surface and in extracellular matrix. The encoded preproprotein is proteolytically processed to generate the mature protein. Fibronectin is involved in cell adhesion and migration processes including embryogenesis, wound healing, blood coagulation, host defense, and metastasis. The gene has three regions subject to alternative splicing, with the potential to produce 20 different transcript variants, at

least one of which encodes an isoform that undergoes proteolytic processing. The full-length

nature of some variants has not been determined. [provided by RefSeq, Jan 2016]

**Protein Families:** Druggable Genome, ES Cell Differentiation/IPS, Secreted Protein

**Protein Pathways:** ECM-receptor interaction, Focal adhesion, Pathways in cancer, Regulation of actin

cytoskeleton, Small cell lung cancer