

## Product datasheet for **TP721182L**

### IL22 (NM\_020525) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human interleukin 22 (IL22)
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Ala34-Ile179
Tag:	Tag Free
Predicted MW:	16.9 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)
Reconstitution Method:	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in ddH <sub>2</sub> O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Storage:	Store at -80°C.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	<a href="#">NP_065386</a>
Locus ID:	50616
UniProt ID:	<a href="#">Q9GZX6</a>
RefSeq Size:	1147
Cytogenetics:	12q15
RefSeq ORF:	537
Synonyms:	IL-21; IL-22; IL-D110; IL-TIF; ILTIF; TIFa; TIFIL-23; zcyto18



[View online »](#)

- Summary:** This gene is a member of the IL10 family of cytokines that mediate cellular inflammatory responses. The encoded protein functions in antimicrobial defense at mucosal surfaces and in tissue repair. This protein also has pro-inflammatory properties and plays a role in in the pathogenesis of several intestinal diseases. [provided by RefSeq, Jul 2018]
- Protein Families:** Druggable Genome, Secreted Protein, Transmembrane
- Protein Pathways:** Cytokine-cytokine receptor interaction, Jak-STAT signaling pathway