

Product datasheet for TP721191

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

YY1 (NM_003403) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Purified recombinant protein of Human YY1 transcription factor (YY1)

Species: Human
Expression Host: E. coli

Expression cDNA Clone

Val221-Gly321

or AA Sequence:

Tag: C-His

Predicted MW: 12.6 kDa

Concentration: lot specific

Purity: >95% as determined by SDS-PAGE and Coomassie blue staining

Buffer: Provided lyophilized from a 0.2 μm filtered solution of 20 mM Tris-HCl, 150 mM NaCl

Endotoxin: Endotoxin level is < 0.1 ng/μg of protein (< 1 EU/μg)

Reconstitution Method: Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the

lyophilized protein in ddH2O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Storage: Store at -80°C.

Stability: Stable for at least 6 months from date of receipt under proper storage and handling

conditions.

RefSeq: NP 003394

 Locus ID:
 7528

 UniProt ID:
 P25490

 RefSeq Size:
 3159

 Cytogenetics:
 14q32.2

 RefSeq ORF:
 1242

Synonyms: DELTA; GADEVS; INO80S; NF-E1; UCRBP; YIN-YANG-1





YY1 (NM_003403) Human Recombinant Protein - TP721191

Summary: YY1 is a ubiquitously distributed transcription factor belonging to the GLI-Kruppel class of zinc

finger proteins. The protein is involved in repressing and activating a diverse number of promoters. YY1 may direct histone deacetylases and histone acetyltransferases to a promoter in order to activate or repress the promoter, thus implicating histone modification

in the function of YY1. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome, Transcription Factors