

Product datasheet for **TP721192L**

FKBP25 (FKBP3) (NM_002013) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human FK506 binding protein 3, 25kDa (FKBP3)
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Met1-Asp224
Tag:	N-His
Predicted MW:	27.3 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)
Storage:	Store at -80°C.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	<u>NP_002004</u>
Locus ID:	2287
UniProt ID:	<u>Q00688</u>
RefSeq Size:	1353
Cytogenetics:	14q21.2
RefSeq ORF:	672
Synonyms:	FKBP-3; FKBP-25; FKBP25; PPlase



[View online »](#)

Summary:

The protein encoded by this gene is a member of the immunophilin protein family, which play a role in immunoregulation and basic cellular processes involving protein folding and trafficking. This encoded protein is a cis-trans prolyl isomerase that binds the immunosuppressants FK506 and rapamycin, as well as histone deacetylases, the transcription factor YY1, casein kinase II, and nucleolin. It has a higher affinity for rapamycin than for FK506 and thus may be an important target molecule for immunosuppression by rapamycin. [provided by RefSeq, Sep 2008]

Protein Families:

Druggable Genome, Stem cell - Pluripotency