

## **Product datasheet for TP721193**

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## ATF1 (NM\_005171) Human Recombinant Protein

## **Product data:**

**Product Type:** Recombinant Proteins

**Description:** Purified recombinant protein of Human activating transcription factor 1 (ATF1)

Species: Human
Expression Host: E. coli

**Expression cDNA Clone** 

or AA Sequence:

Met1-Val271

Tag: C-His

Predicted MW: 30.3 kDa

Concentration: lot specific

**Purity:** >95% as determined by SDS-PAGE and Coomassie blue staining

Buffer: Provided lyophilized from a 0.2 μm filtered solution of 20 mM Tris-HCl, 150 mM NaCl

**Endotoxin:** Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)

**Reconstitution Method:** Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the

lyophilized protein in ddH2O. It is not recommended to reconstitute a concentration less than 100  $\mu$ g/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Storage: Store at -80°C.

Stability: Stable for at least 6 months from date of receipt under proper storage and handling

conditions.

**RefSeq:** NP 005162

Locus ID: 466

UniProt ID: P18846

RefSeq Size: 2505

Cytogenetics: 12q13.12

RefSeq ORF: 813

**Synonyms:** EWS-ATF1; FUS/ATF-1; TREB36



ORÏGENE

**Summary:** 

This gene encodes an activating transcription factor, which belongs to the ATF subfamily and bZIP (basic-region leucine zipper) family. It influences cellular physiologic processes by regulating the expression of downstream target genes, which are related to growth, survival, and other cellular activities. This protein is phosphorylated at serine 63 in its kinase-inducible domain by serine/threonine kinases, cAMP-dependent protein kinase A, calmodulin-dependent protein kinase I/II, mitogen- and stress-activated protein kinase and cyclin-dependent kinase 3 (cdk-3). Its phosphorylation enhances its transactivation and transcriptional activities, and enhances cell transformation. Fusion of this gene and FUS on chromosome 16 or EWSR1 on chromosome 22 induced by translocation generates chimeric proteins in angiomatoid fibrous histiocytoma and clear cell sarcoma. This gene has a pseudogene on chromosome 6. [provided by RefSeq, Aug 2010]

**Protein Families:** 

Druggable Genome, Transcription Factors