

Product datasheet for **TP721206XL**

Hsp22 (HSPB8) (NM_014365) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human heat shock 22kDa protein 8 (HSPB8)
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Met1-Thr196
Tag:	C-His
Predicted MW:	22.7 kDa
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Supplied as a 0.2 um filtered solution of 20mM PB,150mM NaCl,5mM EDTA,pH7.4.
Endotoxin:	Endotoxin level is < 0.1 ng/μg of protein (< 1 EU/μg)
Storage:	Store at < -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_055180
Locus ID:	26353
UniProt ID:	Q9UJY1
RefSeq Size:	2056
Cytogenetics:	12q24.23
RefSeq ORF:	588
Synonyms:	CMT2L; DHMN2; E2IG1; H11; HMN2; HMN2A; HSP22



[View online »](#)

Summary:

The protein encoded by this gene belongs to the superfamily of small heat-shock proteins containing a conservative alpha-crystallin domain at the C-terminal part of the molecule. The expression of this gene is induced by estrogen in estrogen receptor-positive breast cancer cells, and this protein also functions as a chaperone in association with Bag3, a stimulator of macroautophagy. Thus, this gene appears to be involved in regulation of cell proliferation, apoptosis, and carcinogenesis, and mutations in this gene have been associated with different neuromuscular diseases, including Charcot-Marie-Tooth disease. [provided by RefSeq, Jul 2008]

Protein Families:

Druggable Genome, Protein Kinase