

Product datasheet for **TP721213L**

PTPRN (NM_002846) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human protein tyrosine phosphatase, receptor type, N (PTPRN), transcript variant 1
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Gln607-Asn686&Trp795-Leu888
Tag:	N-His
Predicted MW:	22.9 kDa
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_002837
Locus ID:	5798
UniProt ID:	Q16849 , Q96IA0
RefSeq Size:	3649
Cytogenetics:	2q35
RefSeq ORF:	2937
Synonyms:	IA-2; IA-2/PTP; IA2; ICA512; R-PTP-N



[View online »](#)

Summary:

The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitotic cycle, and oncogenic transformation. This PTP possesses an extracellular region, a single transmembrane region, and a single catalytic domain, and thus represents a receptor-type PTP. This PTP was found to be an autoantigen that is reactive with insulin-dependent diabetes mellitus (IDDM) patient sera, and thus may be a potential target of autoimmunity in diabetes mellitus. Alternate splicing results in multiple transcript variants.[provided by RefSeq, Dec 2010]

Protein Families:

Druggable Genome, Transmembrane

Protein Pathways:

Type I diabetes mellitus