

Product datasheet for TP721215L

NRG1 (NM_013962) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins Description: Purified recombinant protein of Human neuregulin 1 (NRG1), transcript variant GGF2 Species: Human E. coli **Expression Host: Expression cDNA Clone** Met1-Lys246 or AA Sequence: N-His Tag: Predicted MW: 29 kDa **Concentration:** lot specific **Purity:** >95% as determined by SDS-PAGE and Coomassie blue staining **Buffer:** Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl Endotoxin: Endotoxin level is < 0.1 ng/ μ g of protein (< 1 EU/ μ g) **Reconstitution Method:** Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Dissolve the lyophilized protein in 4mM HCl. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Store at -80°C. Storage: Stability: Stable for at least 6 months from date of receipt under proper storage and handling conditions. NP 039256 RefSeq: Locus ID: 3084 UniProt ID: Q02297 **RefSeq Size:** 1986 Cytogenetics: 8p12 **RefSeq ORF:** 1266 ARIA; GGF; GGF2; HGL; HRG; HRG1; HRGA; MST131; MSTP131; NDF; NRG1-IT2; SMDF Synonyms:



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

The protein encoded by this gene is a membrane glycoprotein that mediates cell-cell
signaling and plays a critical role in the growth and development of multiple organ systems. An extraordinary variety of different isoforms are produced from this gene through alternative promoter usage and splicing. These isoforms are expressed in a tissue-specific manner and differ significantly in their structure, and are classified as types I, II, III, IV, V and VI. Dysregulation of this gene has been linked to diseases such as cancer, schizophrenia, and bipolar disorder (BPD). [provided by RefSeq, Apr 2016]
Druggable Genome, Secreted Protein, Transcription Factors, Transmembrane
ErbB signaling pathway

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US