

Product datasheet for **TP721227XL**

VSTM2A (NM_182546) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human V-set and transmembrane domain containing 2A (VSTM2A)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Ser24-Phe243
Tag:	C-His
Predicted MW:	25.5 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Endotoxin:	Endotoxin level is < 0.1 ng/μg of protein (< 1 EU/μg)
Reconstitution Method:	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in ddH ₂ O. It is not recommended to reconstitute a concentration less than 100 μg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Storage:	Lyophilized protein should be stored at < -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_872352
Locus ID:	222008
UniProt ID:	Q8TAG5
RefSeq Size:	3122
Cytogenetics:	7p11.2
RefSeq ORF:	729



[View online »](#)

Synonyms: VSTM2

Summary: Plays a role in the regulation of the early stage of white and brown preadipocyte cell differentiation. Promotes adipogenic commitment of preadipocytes by increasing gene expression of the transcription factor PPARG in a BMP4-dependent signaling pathway. [UniProtKB/Swiss-Prot Function]

Protein Families: Secreted Protein, Transmembrane