

## Product datasheet for **TP721253M**

### FCGR1A Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	APC Conjugated Human CD64 Protein (C-His)
Species:	Human
Expression Host:	CHO
Expression cDNA Clone or AA Sequence:	Gln16-Thr287
Tag:	C-His
Predicted MW:	The protein has the predicted molecular weight of 32kDa and migrates at approximately 50-60 kDa on SDS-PAGE with DTT-reduced condition before APC conjugation.
Concentration:	25µg size is bottled at 0.1mg/mL concentration. 100 µg size is bottled at lot specific concentration.
Purity:	>90%
Conjugation:	APC
Buffer:	1xPBS buffer, pH7.4, 0.09% NaN <sub>3</sub> with a carrier protein
Bioactivity:	Positive

The definition of the active protein (purified and biotinylated) is defined as the protein that can bind to its biological receptor/ligand. For conjugated protein, it is defined with its function to bind to the ScFv of the active CAR-transfected cells in flow cytometry test.

Preparation:	Affinity Ni-NTA
Applications:	FACS
Storage:	An unopened vial can be stored at 4°C for 2 weeks or at -20°C and below for six months. This stock solution should be aliquoted and stored at ≤ -70°C to minimize the freeze/thaw cycles.
Stability:	6 Months
RefSeq:	<a href="#">NP_000557.1</a>
Locus ID:	2209
UniProt ID:	<a href="#">P12314</a>



[View online >](#)