

## Product datasheet for **TP721263M**

### CD274 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	APC Conjugated Human PD-L1 Protein (C-Fc)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Phe19-Thr239
Tag:	C-Fc
Predicted MW:	The protein has a predicted molecular weight of 53kDa. The protein migrates at approximately 70kDa on SDS-PAGE with DTT-reduced conditions before APC conjugation.
Concentration:	25µg size is bottled at 0.1mg/mL concentration. 100 µg size is bottled at lot specific concentration.
Purity:	>90%
Conjugation:	APC
Buffer:	1xPBS buffer, pH7.4, 0.09% NaN3 with a carrier protein
Bioactivity:	Positive
	The definition of the active protein (purified and biotinylated) is defined as the protein that can bind to its biological receptor/ligand. For conjugated protein, it is defined with its function to bind to the ScFv of the active CAR-transfected cells in flow cytometry test.
Preparation:	Affinity Protein A
Applications:	FACS
Storage:	An unopened vial can be stored at 4°C for 2 weeks or at -20°C and below for six months. This stock solution should be aliquoted and stored at ≤ -70°C to minimize the freeze/thaw cycles.
Stability:	6 Months
RefSeq:	<a href="#">NP_054862.1</a>
Locus ID:	29126
UniProt ID:	<a href="#">Q9NZQ7</a>



[View online »](#)