

Product datasheet for TP721269

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: APC Conjugated Human Mesothelin (E296-G580) Protein (C-His)

Species:HumanExpression Host:HEK293

Expression cDNA Clone

Expression CDNA Cione

or AA Sequence:

E296-G580

Tag: C-His

Predicted MW: The protein has a predicted molecular weight of 33.7 kDa and migrates at approximately 35-

45 kDa on SDS-PAGE with DTT-reduced conditions before APC conjugation.

Concentration: 25µg size is bottled at 0.1mg/mL concentration. 100 µg size is bottled at lot specific

concentration.

Purity: >90% Conjugation: APC

Buffer: 1xPBS buffer, pH7.4, 0.09% NaN3 with a carrier protein

Bioactivity: Positive

The definition of the active protein (purified and biotinylated) is defined as the protein that can bind to its biological receptor/ligand. For conjugated protein, it is defined with its function

to bind to the ScFv of the active CAR-transfected cells in flow cytometry test.

Preparation: Affinity Ni-NTA

Applications: FACS

Storage: An unopened vial can be stored at 4°C for 2 weeks or at -20°C and below for six months. This

stock solution should be aliquoted and stored at \leq -70°C to minimize the freeze/thaw cycles.

Stability: 6 Months

RefSeq: <u>AAH09272.1</u>

Locus ID: 10232

UniProt ID: Q13421

