

Product datasheet for TP721270

OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

Rockville, MD 20850, US
Phone: +1-888-267-4436
https://www.origene.com
techsupport@origene.com
EU: info-de@origene.com
CN: techsupport@origene.cn

Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Human Mesothelin (E296-G580) Protein (C-Fc)

Species: Human Expression Host: HEK293

Expression cDNA Clone

E296-G580

or AA Sequence:

Tag: C-His

Predicted MW: The protein has a predicted molecular weight of 60.7 kDa and migrates at approximately 70

kDa on SDS-PAGE with DTT-reduced conditions.

Concentration: 25µg size is bottled at 0.2mg/mL concentration. 100 µg size is bottled at lot specific

concentration.

Purity: >90%

Buffer: 1xPBS buffer, pH7.4

Bioactivity: Positive

The definition of the active protein (purified and biotinylated) is defined as the protein that can bind to its biological receptor/ligand. For conjugated protein, it is defined with its function

to bind to the ScFv of the active CAR-transfected cells in flow cytometry test.

Preparation: Affinity Protein A

Applications: ELISA

Storage: An unopened vial can be stored at 4°C for 2 weeks or at -20°C and below for six months. This

stock solution should be aliquoted and stored at \leq -70°C to minimize the freeze/thaw cycles.

Stability: 6 Months

RefSeq: AAH09272.1


Locus ID: 10232

UniProt ID: Q13421





Product images:



Recombinant Human MSLN (E296-G580) Protein (C-Fc) on SDS-PAGE under reducing condition P(+) and non-reducing condition P(-). The purity of this protein appears to be greater than 95% based on Coomassie-blue staining.