

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# Product datasheet for TP721283M

#### **TNFRSF17 Human Recombinant Protein**


### **Product data:**

Product Type:	Recombinant Proteins
Description:	Human BCMA Protein (C-Fc-Avi)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Met1-Ala54
Tag:	C-Fc-Avi
Predicted MW:	The protein has a predicted molecular weight of 34.3 kDa and migrates at approximately 35- 45 kDa on SDS-PAGE with DTT-reduced conditions.
Concentration:	25 $\mu$ g size is bottled at 0.2mg/mL concentration. 100 $\mu$ g size is bottled at lot specific concentration.
Purity:	>90%
Buffer:	1xPBS buffer, pH7.4
Bioactivity:	Positive
	The definition of the active protein (purified and biotinylated) is defined as the protein that can bind to its biological receptor/ligand. For conjugated protein, it is defined with its function to bind to the ScFv of the active CAR-transfected cells in flow cytometry test.
Preparation:	Affinity Protein A
Applications:	ELISA
Storage:	An unopened vial can be stored at 4°C for 2 weeks or at -20°C and below for six months. This stock solution should be aliquoted and stored at ≤ -70°C to minimize the freeze/thaw cycles.
Stability:	6 Months
RefSeq:	<u>NP 001183.2</u>
Locus ID:	608
UniProt ID:	<u>Q02223</u>



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

## **Product images:**



Recombinant Human BCMA Protein (C-Fc-Avi) on SDS-PAGE under reducing condition P(+) and non-reducing condition P(-). The purity of this protein appears to be greater than 95% based on Coomassie-blue staining.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US