

# Product datasheet for TP721290M

## FOLR1 Human Recombinant Protein

### **Product data:**

#### OriGene Technologies, Inc.


9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	Recombinant Proteins
Description:	PE Conjugated Human FOLR1 Protein (C-His)
Species:	Human
Expression Host:	СНО
Expression cDNA Clone or AA Sequence:	Met1-Met233
Tag:	C-His
Predicted MW:	The protein has a predicted molecular weight of 26 kDa and migrates at approximately 35-45 kDa on SDS-PAGE with DTT-reduced conditions before PE conjugation.
Concentration:	25 $\mu$ g size is bottled at 0.1mg/mL concentration. 100 $\mu$ g size is bottled at lot specific concentration.
Purity:	>90%
Conjugation:	PE
Buffer:	1xPBS buffer, pH7.4, 0.09% NaN3 with a carrier protein
Bioactivity:	Positive
	The definition of the active protein (purified and biotinylated) is defined as the protein that can bind to its biological receptor/ligand. For conjugated protein, it is defined with its function to bind to the ScFv of the active CAR-transfected cells in flow cytometry test.
Preparation:	Affinity Ni-NTA
Applications:	FACS
Storage:	An unopened vial can be stored at 4°C for 2 weeks or at -20°C and below for six months. This stock solution should be aliquoted and stored at $\leq$ -70°C to minimize the freeze/thaw cycles.
Stability:	6 Months
RefSeq:	<u>NP 057937.1</u>
Locus ID:	2348
UniProt ID:	<u>P15328</u>



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US 

### **Product images:**



PE conjugated Human FOLR1 Protein (C-His) on SDS-PAGE under reducing condition PE (+) and non-reducing condition PE(-). C(-) is unconjugated Human FOLR1 Protein under non reducing condition.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US