

## **Product datasheet for TP721293M**

# **OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

### **Human Recombinant Protein**

**Product data:** 

**Product Type:** Recombinant Proteins

**Description:** Human PD1 Protein (C-His-Avi)

Species: Human Expression Host: HEK293

**Expression cDNA Clone** 

Leu25-Gln167

or AA Sequence:

\_\_\_\_\_\_

Tag: C-His-Avi

Predicted MW: The protein has a predicted molecular weight of 28 kDa and migrates at approximately 30-45

kDa on SDS-PAGE with DTT-reduced conditions.

**Concentration:** 25μg size is bottled at 0.2mg/mL concentration. 100 μg size is bottled at lot specific

concentration.

**Purity:** >90%

**Buffer:** 1xPBS buffer, pH7.4

**Bioactivity:** Positive

The definition of the active protein (purified and biotinylated) is defined as the protein that can bind to its biological receptor/ligand. For conjugated protein, it is defined with its function

to bind to the ScFv of the active CAR-transfected cells in flow cytometry test.

**Preparation:** Affinity Ni-NTA

**Applications:** ELISA

**Storage:** An unopened vial can be stored at 4°C for 2 weeks or at -20°C and below for six months. This

stock solution should be aliquoted and stored at  $\leq$  -70°C to minimize the freeze/thaw cycles.

**Stability:** 6 Months

**RefSeq:** NP 005009.2


**Locus ID:** 5133

UniProt ID: Q15116





# **Product images:**



Human PD1 Protein (C-His-Avi) on SDS-PAGE under reducing condition P(+) and non-reducing condition P(-). The purity of this protein appears to be greater than 95% based on Coomassie-blue staining.