

Product datasheet for **TP721309M**


CD8A Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Human CD8 α Protein (C-His-Avi)
Species:	Human
Expression Host:	CHO
Expression cDNA Clone or AA Sequence:	Ser22-Ala182
Tag:	C-His-Avi
Predicted MW:	The protein has a predicted molecular weight of 21 kDa and migrates at approximately 25-30 kDa on SDS-PAGE with DTT-reduced conditions.
Concentration:	25 μ g size is bottled at 0.2mg/mL concentration. 100 μ g size is bottled at lot specific concentration.
Purity:	>90%
Buffer:	1xPBS buffer, pH7.4
Bioactivity:	Positive
	<p>The definition of the active protein (purified and biotinylated) is defined as the protein that can bind to its biological receptor/ligand. For conjugated protein, it is defined with its function to bind to the ScFv of the active CAR-transfected cells in flow cytometry test.</p>
Preparation:	Affinity Ni-NTA
Applications:	ELISA
Storage:	An unopened vial can be stored at 4°C for 2 weeks or at -20°C and below for six months. This stock solution should be aliquoted and stored at \leq -70°C to minimize the freeze/thaw cycles.
Stability:	6 Months
RefSeq:	NP_001759.3
Locus ID:	925
UniProt ID:	P01732



[View online »](#)

Product images:

Human CD8α Protein (C-His-Avi) on SDS-PAGE under reducing condition P(+) and non-reducing condition P(-). The purity of this protein appears to be greater than 95% based on Coomassie-blue staining.