

Product datasheet for TP721316

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436

OriGene Technologies, Inc.

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

ERBB2 Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Biotinylated Human HER2 Protein (C-His-Avi)

Species: Human Expression Host: CHO

Expression cDNA Clone

Ser22-Thr652

or AA Sequence:

Tag:

C-His-Avi

Predicted MW: The protein has a predicted molecular weight of 73 kDa and migrates at approximately 80

kDa on SDS-PAGE with DTT-reduced conditions.

Concentration: 25µg size is bottled at 0.2mg/mL concentration. 100 µg size is bottled at lot specific

concentration.

Purity: >90%

Conjugation: Biotin

Buffer: 1xPBS buffer, pH7.4

Bioactivity: Positive

The definition of the active protein (purified and biotinylated) is defined as the protein that can bind to its biological receptor/ligand. For conjugated protein, it is defined with its function

to bind to the ScFv of the active CAR-transfected cells in flow cytometry test.

Preparation: Affinity Ni-NTA

Applications: ELISA

Storage: An unopened vial can be stored at 4°C for 2 weeks or at -20°C and below for six months. This

stock solution should be aliquoted and stored at \leq -70°C to minimize the freeze/thaw cycles.

Stability: 6 Months

RefSeq: <u>NP 004439.2</u>

Locus ID: 2064 UniProt ID: P04626






Summary:

Human epithelial growth factor receptor 2 (HER2), also named ErbB2/Neu receptor, is a member of the epidermal growth factor receptor (EGFR; also known as ErbB) family of receptor tyrosine kinases. In human, this family are consisted of four members: HER1 (EGFR, ERBB1), HER2 (ERBB2), HER3 (ERBB3) and HER4 (ERBB4). The HER family proteins are type I transmembrane growth factor receptors that function to activate intracellular signaling pathways in response to extracellular signals. Their structure consists of an extracellular ligand binding domain, a transmembrane domain, and an intracellular tyrosine kinase domain. Unlike other members of the family, HER2 lacks ligand binding activity and its signaling function is engaged by its ligand-bound heterodimeric partners. Its expression has a close relationship with various tumors. Its overexpression is found in malignant tumors, such as breast, ovarian, gastric, and colorectal cancers.

Product images:



Biotinylated Human HER2 Protein (C-His-Avi) on SDS-PAGE under reducing condition P(+) and non-reducing condition P(-). The purity of this protein appears to be greater than 90% based on Coomassie-blue staining.