

Product datasheet for TP721329M

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com

OriGene Technologies, Inc.

EU: info-de@origene.com CN: techsupport@origene.cn

L1CAM Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Human CD171/L1CAM Protein (C-His)

Species: Human Expression Host: CHO

Expression cDNA Clone

lle20-Glu1120

or AA Sequence:

Tag:

C-His

Predicted MW: The protein has a predicted molecular weight of 125 kDa and migrates at approximately 160-

180 kDa on SDS-PAGE with DTT-reduced conditions.

Concentration: 25µg size is bottled at 0.2mg/mL concentration. 100 µg size is bottled at lot specific

concentration.

Purity: >90%

Buffer: 1xPBS buffer, pH7.4

Bioactivity: Positive

The definition of the active protein (purified and biotinylated) is defined as the protein that can bind to its biological receptor/ligand. For conjugated protein, it is defined with its function

to bind to the ScFv of the active CAR-transfected cells in flow cytometry test.

Preparation: Affinity Ni-NTA

Applications: ELISA

Storage: An unopened vial can be stored at 4°C for 2 weeks or at -20°C and below for six months. This

stock solution should be aliquoted and stored at \leq -70°C to minimize the freeze/thaw cycles.

Stability: 6 Months

RefSeq: <u>NP 000416.1</u>


Locus ID: 3897

UniProt ID: P32004





Product images:



Human CD171 Protein (C-His) on SDS-PAGE under reducing condition P(+) and non-reducing condition P(-). The purity of this protein appears to be greater than 90% based on Coomassie-blue staining.