

Product datasheet for TP721339

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436

OriGene Technologies, Inc.

Rockville, MD 20850, US
Phone: +1-888-267-4436
https://www.origene.com
techsupport@origene.com
EU: info-de@origene.com
CN: techsupport@origene.cn

Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Human Trop2 Protein (C-His)

Species: Human Expression Host: CHO

Expression cDNA Clone

Gln31-Thr274

or AA Sequence:

Tag: C-His

Predicted MW: The protein has a predicted molecular weight of 31kDa and migrates at approximately 40-

50kDa on SDS-PAGE with DTT-reduced condition.

Concentration: 25µg size is bottled at 0.2mg/mL concentration. 100 µg size is bottled at lot specific

concentration.

Purity: >90%

Buffer: 1xPBS buffer, pH7.4

Bioactivity: Positive

The definition of the active protein (purified and biotinylated) is defined as the protein that can bind to its biological receptor/ligand. For conjugated protein, it is defined with its function

to bind to the ScFv of the active CAR-transfected cells in flow cytometry test.

Preparation: Affinity Ni-NTA

Applications: ELISA

Storage: An unopened vial can be stored at 4°C for 2 weeks or at -20°C and below for six months. This

stock solution should be aliquoted and stored at \leq -70°C to minimize the freeze/thaw cycles.

Stability: 6 Months

RefSeg: NP 002344.2

Locus ID: 4070

UniProt ID: P09758






Summary:


Human TROP-2 protein is a type I cell surface glycoprotein highly expressed in a wide variety of human carcinomas. It is also referred to as tumor-associated calcium signal transducer 2 (TACSTD2), gp50, GA733-1 or M1S1. TROP-2 transduces an intracellular calcium signal, which are consistent with the hypothesis that it acts as a cell surface receptor and support a search for a physiological ligand.

Product images:



Human Trop2 protein (C-His) is coated at 100ng/well. Anti-Trop2 antibody can detect the Trop2 protein with a linear range of 0.02 -2.5 ng/mL

TROP-2 ELISA



Human Trop2 Protein (C-His) on SDS-PAGE under reducing condition P(+) and non-reducing condition P(-). The purity of this protein appears to be greater than 95% based on Coomassie-blue staining.