

Product datasheet for TP721342M

Human Recombinant Protein

Product data:

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	Recombinant Proteins
Description:	PE Conjugated Human Trop2 Protein (C-His)
Species:	Human
Expression Host:	СНО
Expression cDNA Clone or AA Sequence:	Gln31-Thr274
Tag:	C-His
Predicted MW:	The protein has a predicted molecular weight of 31kDa and migrates at approximately 40- 50kDa on SDS-PAGE with DTT-reduced condition before PE conjugation.
Concentration:	25 μg size is bottled at 0.1mg/mL concentration. 100 μg size is bottled at lot specific concentration.
Purity:	>90%
Conjugation:	PE
Buffer:	1xPBS buffer, pH7.4, 0.09% NaN3 with a carrier protein
Bioactivity:	Positive
	The definition of the active protein (purified and biotinylated) is defined as the protein that can bind to its biological receptor/ligand. For conjugated protein, it is defined with its function to bind to the ScFv of the active CAR-transfected cells in flow cytometry test.
Preparation:	Affinity Ni-NTA
Applications:	FACS
Storage:	An unopened vial can be stored at 4°C for 2 weeks or at -20°C and below for six months. This stock solution should be aliquoted and stored at \leq -70°C to minimize the freeze/thaw cycles.
Stability:	6 Months
RefSeq:	<u>NP 002344.2</u>
Locus ID:	4070
UniProt ID:	<u>P09758</u>



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US