

Product datasheet for **TP721374M**

LDHA Mouse Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant mouse Lactate Dehydrogenase A/LDHA protein
Species:	Mouse
Expression Host:	Baculovirus
Expression cDNA Clone or AA Sequence:	1-332 aa
Tag:	His-Tag
Predicted MW:	37.5 kDa (340aa)
Concentration:	0.5mg/ml (determined by absorbance at 280nm)
Purity:	> 95% by SDS-PAGE
Buffer:	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol
Storage:	Can be stored at +2°C to +8°C for 1 week. For long term storage, aliquot and store at -20°C to -80°C. Avoid repeated freezing and thawing cycles.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_034829
Locus ID:	3939
UniProt ID:	Q564E2 , P06151



[View online »](#)

Product images:



Coomassie Blue staining of purified recombinant mouse Lactate Dehydrogenase A/LDHA protein. (Cat# [TP721374]).