

Product datasheet for **TP721376M**

IL17a Mouse Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant mouse IL-17A protein
Species:	Mouse
Expression Host:	Baculovirus
Expression cDNA Clone or AA Sequence:	26-158 aa
Tag:	His-Tag
Predicted MW:	16kDa (141aa)
Concentration:	0.5mg/ml (determined by absorbance at 280nm)
Purity:	> 95% by SDS-PAGE
Buffer:	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol
Storage:	Can be stored at +2°C to +8°C for 1 week. For long term storage, aliquot and store at -20°C to -80°C. Avoid repeated freezing and thawing cycles.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_034682.1
Locus ID:	3605
UniProt ID:	Q16552




[View online »](#)

Product images:



Coomassie Blue staining of purified recombinant mouse IL-17A protein. (Cat# [TP721376]).



Recombinant mouse IL-17A protein stimulates IL-6 secretion of NIH3T3 cells. The ED50 range ≤ 2 ng/ml.