

Product datasheet for **TP721377L**


HAPLN1 Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant human HAPLN1 protein
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	16-354 aa
Tag:	His-Tag
Predicted MW:	39.3kDa (345aa)
Concentration:	0.25mg/ml (determined by Absorbance at 280nm)
Purity:	> 95% by SDS-PAGE
Buffer:	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol
Storage:	Can be stored at +2°C to +8°C for 1 week. For long term storage, aliquot and store at -20°C to -80°C. Avoid repeated freezing and thawing cycles.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001875.1
Locus ID:	1404
UniProt ID:	P10915



[View online »](#)

Product images:


Coomassie Blue staining of purified recombinant human HAPLN1 protein. (Cat# [TP721377]).



Measured by binding ability of recombinant human HAPLN1 protein in a functional ELISA with Hyaluronic acid. The ED50 range \leq 1 ug/ml.