

Product datasheet for **TP721845**

TA2R9 Human Recombinant Protein, Synthetic Nanodisc

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human TA2R9 full length protein-Synthetic Nanodisc, 10µg
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	C-terminal Flag tagged overexpression cDNA clone
Tag:	C-term Flag Tag
Predicted MW:	The human full length TA2R9 protein has a MW of 35.6kDa
Concentration:	Please refer to the Certificate of Analysis (COA) for the lot-specific concentration before lyophilization.
Purity:	The purity of the protein is determined by SDS-PAGE and Coomassie blue staining.
Buffer:	Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Do not use solvents with pH lower than 6.5 in subsequent experiments.
Preparation:	Affinity purification
Storage:	Store at -20°C to -80°C for 12 months in lyophilized form.
Stability:	After reconstitution, if not intended for use within a month, aliquot and store at -80°C . Avoid repeated freezing and thawing.
Locus ID:	50835
UniProt ID:	Q9NYW1
Synonyms:	T2R9, TRB6
Protein Families:	ES Cell Differentiation/IPS ,Transmembrane,Stem cell relevant signaling - DSL/Notch pathway,Druggable Genome
Protein Pathways:	GPCRDB Other,Hedgehog Netpath 10,Hedgehog Netpath 10,Apoptosis,Cancer,Notch



[View online »](#)