

Product datasheet for **TP723022**

Enzymogenes Arg-C (Ier) Lysobacter Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Lysobacter Enzymogenes Arg-C (Ier)
Species:	Lysobacter
Expression Host:	Hi-5 insect
Tag:	C-His
Predicted MW:	26.8 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Lyophilized from a 0.2 μ M filtered solution of 20mM phosphate buffer, 100mM NaCl, pH 7.2
Bioactivity:	The reaction is measured as an increase in absorbance at 253 nm resulting from the hydrolysis of N-benzoyl-L-arginine ethyl ester (BAEE).
Endotoxin:	Endotoxin level is < 0.1 ng/ μ g of protein (< 1 EU/ μ g)
Storage:	Store at -80°C.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	AAD11571



[View online »](#)