


Product datasheet for TP723023

ADI Mycoplasma Recombinant Protein (ADI) Mycoplasma Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mycoplasma Arginine Deiminase (ADI)
Species:	Mycoplasma
Expression Host:	E. coli
Tag:	Tag Free
Predicted MW:	46.3 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Lyophilized from a 0.2 μ M filtered solution of 20mM phosphate buffer, 100mM NaCl, pH 7.2
Bioactivity:	Measured by its ability to induce apoptosis in Jurkat cells using a concentration of 100-150 ng/ml.
Endotoxin:	Endotoxin level is < 0.1 ng/ μ g of protein (< 1 EU/ μ g)
Storage:	Store at -80°C.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	CAA38210

Product images:



[View online »](#)