

## Product datasheet for **TP723141**

### Glu-C *Escherichia coli* Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Staphylococcus Glu-C
Species:	Escherichia coli
Expression Host:	E. coli
Tag:	Tag Free
Predicted MW:	14.6 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Lyophilized from a 0.2 $\mu$ M filtered solution of 20mM phosphate buffer, 100mM NaCl, pH 7.2
Bioactivity:	Cleaves at the Carboxyl side of E (can also cleave D under certain conditions).
Endotoxin:	Endotoxin level is < 0.1 ng/ $\mu$ g of protein (< 1 EU/ $\mu$ g)
Storage:	Store at -80°C.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	<a href="#">CAA68434</a>

#### Product images:

