

<b>Name: Purified recombinant protein of Human neuregulin 1 (NRG1), transcript variant GGF2.</b> <b>Product Data Sheet - PROTEIN</b>	<b>Catalog: TP723155</b>
---	--------------------------

**Gene Name:** neuregulin 1  
**GeneBank accession:** NM\_013962

**Description:** Purified recombinant protein of Human neuregulin 1 (NRG1), transcript variant GGF2. (TP723155)

**Amount:** 50 ug

**Buffer and Storage:** Lyophilized from a 0.2 μM filtered solution of 20mM phosphate buffer,100mM NaCl, pH 7.2

**Source:** Recombinant protein was produced with  
 SHLVKCAEKEKTFVNGGECFMVKDLSNPSRYLCKCPNEFTGDRCQNYVMASFYKHLGIE  
 FMEAE

**Protein Accn:** NP\_001153467

**Gene Synonym:** ARIA|GGF|GGF2|HGL|HRG|HRG1|HRGA|MST131|MSTP131|NDF|NRG1-IT2|S  
 MDF


**Predicted molecular weight:** 7.5 kDa

**Tags:** Tag-free

**Purity:** >95% as determined by SDS-PAGE and Coomassie blue staining

**Concentration:** Resuspend the protein in the desired concentration in proper buffer

**Protein Gel Analysis Data:**



Purified recombinant protein NRG1 was analyzed by SDS-PAGE gel and Coomassie Blue Staining.

\* More data may be available on our website: [TP723155](#)