


Name: Purified recombinant protein of Rat interleukin 4 (Il4). Product Data Sheet - PROTEIN	Catalog: TP723236
--	--------------------------

Gene Name: interleukin 4
GeneBank accession: NM_201270

Description: Purified recombinant protein of Rat interleukin 4 (Il4). (TP723236)
Amount: 20 ug
Buffer and Storage: Lyophilized from a 0.2 μM filtered solution of 20mM phosphate buffer,100mM NaCl, pH 7.2
Source: Recombinant protein was produced with
 MHGCNDSPLREIINTLNQVTEKGPCTEMFVPDVLATRNTTENELICRASRVLRKFYFPR
 DVPPCLKNKSGVLGELRKLCRGVSGLNLSRCTVNESTLTTLKDFLESLSILRGKYLQSCTS
 MS
Protein Accn: NP_958427
Gene Synonym: Il4e12
Predicted molecular weight: 14 kDa
Tags: Tag-free
Purity: >95% as determined by SDS-PAGE and Coomassie blue staining
Concentration: Resuspend the protein in the desired concentration in proper buffer

Protein Gel Analysis Data:



Purified recombinant protein Il4 was analyzed by SDS-PAGE gel and Coomassie Blue Staining.

* More data may be available on our website: [TP723236](https://www.origene.com/TP723236)